Main Features



1. Prevention of Patient (Paed) Abduction


- By alerting the staffs on the dashboard when...
 - Safety band has been tampered/removed without authorization or is Undetected
 - Patient Leaves a designated area or detected at the restricted exit without authorization
 - Patient does not return after an authorized period is over
- Auto door lock when patient is near by the door if no authorization give.

Patient Movement Monitoring

 Show the patient latest movement on dashboard when patient is moving within the RFID covered area

3. Patient Body Temperature Monitoring

 Show the patient latest body temperature reading on dashboard when patient is attached with Wireless Temperature tag Smartsense Patient Tracking and Wireless Temperature Dashboard at Paediatric Ward



Two Types of Alerts



1. Physiological (Temperature) Alert 🀠




High Temperature Alert

WD5555-5P2655-19
35.8

Ward 55
TEST FOR XY FOR SAP


Low Temperature Alert

2. Patient Alert




Ward 55

TEST 333





Strap Tampered Alert



Low Battery Alert

Sensor Fault Alert